



## **ALLEGRO GENERAL PURPOSE WITH HEAD SUPPORT SLING**

**\$335.00 - \$459.00** Price range: \$335.00 through \$459.00



**SKU:** GPHS0-H

## **PRODUCT DESCRIPTION**

### **PURPOSE & BENEFITS**

Allegro General Purpose Slings are designed to be simple and safe to use for both home and institutional patient lifting. The slings are manufactured from polyester fabric or mesh (suitable for bathing) with internal padding in pressure areas and utilise webbing attachment loops. The webbing loops are colour coded for sizing to allow easy sling selection and identification. Webbing loops of varying lengths are provided at each connection point to allow a lay back posture for the patient or, by using the shortest loops, a higher lift.

### **FEATURES**

- Sling is easily fitted to a patient in bed support for patients with some or a chair, providing excellent upper upper body control. body and head support.
- Recommended for lifting from the floor and general lifting of the less secure patient.
- Safe Working Load: 300kg

### **SIZING**