



ALLEGRO REFURBISHED SIT TO STAND AID



SKU: RF SAL01

PRODUCT DESCRIPTION

Purpose & Benefits

The **Allegro Refurbished Sit to Stand Aid** is a reliable patient transfer solution designed to assist individuals with limited mobility in safely transitioning from a seated position to a supported standing position. Ideal for home care, aged care facilities, hospitals, and rehabilitation environments, this powered standing transfer aid helps reduce manual handling while promoting safer and more comfortable patient transfers.

Professionally refurbished by Medimart, each unit is thoroughly inspected, cleaned, sanitised, tested, serviced where required, and confirmed to be in full working condition.

Featuring an ergonomic design for both patients and caregivers, the Allegro Sit to Stand Aid provides stable support throughout the transfer process. The electric leg spread allows the base to be easily adjusted for improved access around chairs, wheelchairs, and beds, while the height-adjustable knee pad assists with correct patient positioning. A removable footplate offers additional flexibility to suit different transfer requirements.

Powered by a 24 V DC lead acid battery and capable of supporting users up to **200 kg**, the Allegro Sit to Stand Aid delivers dependable performance for everyday patient handling.

KEY FEATURES

- Professionally refurbished by Medimart
- Designed for safe and controlled sit-to-stand transfers
- Suitable for home care, hospitals, rehabilitation and aged care
- Electric leg spread for easy base adjustment
- Ergonomic design for patient comfort and caregiver usability
- Height-adjustable knee pad for improved positioning
- Removable footplate for greater transfer flexibility
- Battery-powered operation
- Stable and easy to manoeuvre
- Safe Working Load: **200 kg**

SPECIFICATIONS

Specification	Value
Safe Working Load	200 kg
Base Length	101 cm
External Width (Closed)	66.5 cm
External Width (Open)	103 cm
Internal Width (Closed)	56.5 cm
Internal Width (Open)	93 cm
Under Bed Clearance	13 cm
Minimum Lift Arm Height	88 cm
Maximum Lift Arm Height	168.5 cm
Under Base Height	5.5 cm
Lifter Weight	51 kg
Minimum Storage Height	119 cm
Battery	24 V DC, 2.9 Ah Lead Acid

REFURBISHED CONDITION

This product has been professionally refurbished by **Medimart**.

Each unit has been:

- Fully inspected, tested, and **Test & Tagged** to ensure safe operation
- Cleaned and sanitised

- Serviced where required
- Confirmed to be in full working condition

As a refurbished product, the item may display cosmetic signs of previous use, including minor scratches, scuff marks, paint chips, or other surface imperfections. These are normal for refurbished equipment and do not affect the product's safety, functionality, or performance.

The condition of each unit is reflected in the selling price, providing excellent value compared to purchasing a new product.

WARRANTY

This refurbished product is supplied with a **3-month Medimart warranty**.

Terms & Conditions: Medimart Refurbished Product Terms & Conditions apply.

Please Note: Product images are provided as a guide to the general appearance of the refurbished item. While we endeavour to accurately represent the product, images may not show every cosmetic mark or surface imperfection. Cosmetic signs of previous use may vary slightly between units. All refurbished products are professionally cleaned, sanitised, inspected, tested, and confirmed to be in full working condition before sale.

***Please note: Refurbished items are not eligible for free delivery.**

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au