



ALLEGRO STANDARD YOKE WITH INTEGRATED WEIGH SCALE SPREADER BAR

\$2,606.00



SKU: SY02WD

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Allegro Concepts weigh scale option is a universal weighing module that is useful for patients requiring constant monitoring of weight for general wellbeing:

- 320kg SWL
- 100g graduations
- Tare function
- Kg / Lb conversion
- Suitable for yoke or pivot attachments
- Useful where weighing of patients is difficult due to the level of support required

SPECIFICATIONS

- **Compatible with:** [Alto Lifter](#), [Sonata](#) & [Sonata ELS](#)
- All day-to-day patient lifting
- Monitor patient wellbeing
- Minimal loss of hoist lifting height (approx. 20mm)
- 320 kg SWL
- Easily adaptable to any style of patient hoist
- Use with general purpose (loop attachment) slings
- Non-absorbent polyurethane padding