


ALLEGRO TRI HOOK YOKE SPREADER BAR



SKU: SY03



PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Allegro Concepts Hoist Accessories are designed and manufactured in Australia, with a view to providing healthcare facilities with customised patient transfer packages that meet the needs of each individual clients.

- 6 Sling attachment points at extremity of curved yoke bar
- Use with general purpose (loop attachment) slings
- Non-absorbent polyurethane padding

SPECIFICATIONS

- **Compatible with:** [Alto Lifter](#), [Sonata](#) & [Sonata ELS](#)
- **Indications:** Where the sling prescribed has more than the usual 4 attachment loops
- **Useful for:** Facilities which have historically trained with this configuration

For further information go to **[medimart.com.au](https://www.medimart.com.au)** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North

Brighton SA 5048

08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice

Ave, Elizabeth Vale SA 5112

08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,

Reynella SA 5161

08 7079 8771

reynella@medimart.com.au

VIC**Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT**Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au