




## **ALLEGRO WALKING SLING**

**\$376.00**



**SKU: WS0-H**

## **PRODUCT DESCRIPTION**

### **PURPOSE & BENEFITS**

Walking Slings help to provide support and assist when walking in the rehabilitation environment.

Used with ceiling track hoists and compatible floor hoists, these slings provide a safe and supportive aid to assist with those with walking difficulties.

Designed to be simple and safe to use for both home and institutional patient use.

Manufactured in fully quilted fabric, preventing fabric buckling making this sling comfortable and easy to use.

### **FEATURES**

- This sling helps support and assist when walking in the rehabilitation environment.
- Used with Ceiling Track Hoists and Compatible Standing Hoists.
- Safe Working Load: 200kg

### **SIZING**