

## ALRICK 2400 HOME CARE BED BARIATRIC



**SKU:** 7-24HC-D-01101400-IS



## PRODUCT DESCRIPTION

### Purpose & Benefits

The **Alrick 2400 Home Care Bariatric Bed** is an Australian-made electric bed designed to provide advanced comfort, safety and support for bariatric patients in home or care environments. With a Safe Working Load of up to 400 kg, it offers superior functionality with smooth vertical lift, auto-chair positioning and under-bed lighting, while maintaining a modern, stylish appearance suited to residential care.

### Features

- Australian made and owned for local quality assurance and support
- Vertical lift system for improved patient transfer and caregiver ergonomics
- Electric auto-regression backrest and electric knee break for smooth adjustment
- Manual vascular support to assist circulation and comfort
- Forward and reverse Trendelenburg functions for flexible positioning
- Auto chair positioning for user comfort and carer convenience
- Safe Working Load between 350–400 kg for bariatric applications

- 125 mm twin-wheel braking castors for stability and mobility
- Pull-out extension to accommodate taller users
- Back-lit handset for night-time operation and visibility
- Bluetooth connectivity (handset sold separately)
- Battery backup ensures continued function during power outages
- Under-bed lighting for user safety during night hours
- Curved Natural Oak end boards paired with an Ironstone powder-coated frame for a modern finish

## **Optional Features**

- CPR backrest release for emergency situations
- Electric vascular support for enhanced comfort and circulation
- Central locking castor system for added stability
- Bluetooth handset compatibility
- Attendant control for carer-assisted adjustments
- USB charging ports for convenience
- Custom powder coat colour options
- Custom end board style and finish selections
- Extended sleeping surface length (2 m+)

## **Specifications**

**Specification**
**Details**

|                         |                                                                                  |
|-------------------------|----------------------------------------------------------------------------------|
| Model                   | Alrick 2400 Home Care Bariatric Bed                                              |
| Safe Working Load (SWL) | 350 – 400 kg                                                                     |
| Lift Type               | Vertical lift system                                                             |
| Functions               | High/Low, Backrest, Knee Break, Trendelenburg, Reverse Trendelenburg, Auto Chair |
| Castors                 | 125 mm twin-wheel braking castors                                                |
| Frame Finish            | Ironstone powder-coated frame                                                    |
| End Boards              | Curved Natural Oak finish                                                        |
| Controls                | Back-lit handset (Bluetooth compatible, sold separately)                         |
| Additional Features     | Under-bed light, battery backup, pull-out bed extension                          |
| Made In                 | Australia                                                                        |