




Alrick

ALRICK 3004 SERIES PAEDIATRIC CHANGE TABLE



SKU: 7-3004

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Alrick 3004 Series Single Tower Paediatric Change table is ideal for children and small adults up to 90kg. The foot control operation provides flexible and easy-to-control height adjustment, while the single tower lifting column conceals electrics within the column for safe and smooth hi-low functioning.

Moreover, its durable design ensures longevity, making it an economical choice for hospitals and care facilities seeking an effective solution for paediatric patient handling. With its advanced safety features and ergonomic design, the Alrick 3004 Series Single Tower Paediatric Change table promotes a comfortable, secure environment that's conducive to child-centred care.

FEATURES

- Electric Hi-Low
- Single Column
- Foot Control Operation
- Standard 125mm Braking Castors
- Charcoal Padded Vinyl Top
- APO/Ironstone Powder Coated Frame

SPECIFICATIONS

- Overall Dimensions: 1400mm x 700mm
- Height Range: 480mm - 880mm
- Weight Capacity: 90kg