



Alrick

ALRICK 3004 SERIES PAEDIATRIC CHANGE TABLE



SKU: 7-3004

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Alrick 3004 Series Single Tower Paediatric Change table is ideal for children and small adults up to 90kg. The foot control operation provides flexible and easy-to-control height adjustment, while the single tower lifting column conceals electrics within the column for safe and smooth hi-low functioning.

Moreover, its durable design ensures longevity, making it an economical choice for hospitals and care facilities seeking an effective solution for paediatric patient handling. With its advanced safety features and ergonomic design, the Alrick 3004 Series Single Tower Paediatric Change table promotes a comfortable, secure environment that's conducive to child-centred care.

FEATURES

- Electric Hi-Low
- Single Column
- Foot Control Operation
- Standard 125mm Braking Castors
- Charcoal Padded Vinyl Top
- APO/Ironstone Powder Coated Frame

SPECIFICATIONS

- Overall Dimensions: 1400mm x 700mm
- Height Range: 480mm - 880mm
- Weight Capacity: 90kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au