



ALRICK EN9 FLOORLINE BED



SKU: 7-EN9

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The **Alrick EN9 Floorline Bed** is a versatile, Australian-made universal care bed designed for low, high, and bariatric care environments. With an ultra-low height of just 100 mm, it enhances safety for dementia and elderly patients by reducing fall risks while maintaining full functionality for acute and bariatric settings. This three-in-one model features advanced electric adjustments, intelligent ergonomic design, and customisable options for both patient comfort and caregiver efficiency.

FEATURES

- Australian made and owned
- Floorline design with vertical lift for enhanced safety and convenience
- Electric auto-regression backrest and knee-break functions
- Electric and manual vascular support options
- CPR backrest release for emergency use
- Dual locking castors for secure positioning
- Bluetooth handset (HC40) and backlit attendant control
- Battery backup ensures continuous operation during power outages
- Underbed lighting with warm white or custom colour options
- USB charging ports for patient convenience
- Custom powdercoat and endboard colour options
- Luxa full surround available in timber or upholstered finish
- Compatible with floorline partner bed setup
- Ironstone powder-coated frame for durability and style
- Custom sleeping surface length (2m+ available)

SPECIFICATIONS

Size	External Dimensions	Sleeping Surface	Hi-Lo Position	Backrest / Knee Angle	TR / ATR Angle	SWL
Single	920 × 2115 mm	900 × 2000 mm	100 – 725 mm	62° / 35°	15° / 15°	200 kg
King Single	1090 × 2115 mm	1070 × 2000 mm	100 – 725 mm	62° / 35°	15° / 15°	200 kg
Extra Wide	1220 × 2115 mm	1200 × 2000 mm	170 – 725 mm	62° / 35°	15° / 15°	250 kg
Double	1390 × 2115 mm	1370 × 2000 mm	170 – 725 mm	62° / 35°	15° / 15°	250 kg
Queen Double	1390 × 2115 mm	1370 × 2000 mm	170 – 725 mm	62° / 35°	15° / 15°	250 kg
King Double	1820 × 2115 mm	1800 × 2000 mm	170 – 725 mm	62° / 35°	15° / 15°	250 kg

[Click to view the brochure](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au