



AMBASSADOR PREMIER A3 DUAL MOTOR LIFT RECLINER CHAIR



SKU: PA3-H

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

For the most detailed look at the Ambassador recliner please click below to view our brochure.

[Click here to view our brochure.](#)

Ambassador is a market leader in the manufacture and supply of superior lift & recline chairs for the healthcare market.

A great blend of exceptional comfort and advanced technologies, designing the ultimate bariatric recliner chair with exclusive mechanisms and components that are extra strong and durable to provide greater lifting capacity than any other lift and recliner chair in this category.

Ambassador chairs enable more independence in the home and numerous benefits for the caregiver or nursing staff.

FEATURES

(Across all Ambassador bariatric recliner chair range)

- Sit, recline, rest with elevation and stand at the simple touch of a button
- Rechargeable back up battery operating system
- Simultaneous movement of backrest and footrest(A2 and A3 models)
- Elevated footrest above seat height, varying foot and backrest options – ideal for circulation and recovery
- Removable chair covers, seat and backrest cushions for ease of cleaning or a whole new look
- A combination of quality fabrics to personalise the chair to your individual taste
- A selection of width, height, depth and lumbar pillows for the perfect fit

Storage pockets

The A3 – Key Features: Independent footrest & backrest. Leg rest 85 mm higher than the body. Heavy duty TRIPLE motor design. Extra large size. 2 cushions.

Also available Premier A1 (single motor action) with weight capacity 140 kg.

SPECIFICATIONS

Overall width:	1200mm
Overall height:	1150mm
Overall depth:	1000mm
Seat height:	500mm
Seat depth:	520mm
Seat width:	690mm
Armrest height:	190mm
Backrest height:	720mm
Max weight capacity:	250kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriiffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au