


AMBASSADOR PREMIER A4 BARIATRIC RECLINER

\$7,029.75



SKU: PA4MBB-H



PRODUCT DESCRIPTION

Purpose and Benefits

The Ambassador Premier A4 Four Motor Lift and Recline Chair is a high-comfort, high-function seating solution designed for individuals with limited mobility. With a safe working load of 370 kg, this chair supports users with generous stability while assisting with safe transitions between sitting and standing.

This version features a fixed seat width of 69 cm and offers exceptional comfort and durability for daily use. With smooth motorised positioning and user-friendly controls, it is ideal for home environments or residential care where reliability, safety, and ease of use are essential.

Features

- Four-motor system provides independent operation of footrest and backrest
- Elevating leg rest rises 85 mm above seat height to encourage circulation
- Heavy-duty bolted steel frame for dependable long-term support
- Fixed seat width of 69 cm suits most users comfortably

- Soft-touch, breathable fabric upholstery with removable covers for easy maintenance
- Built-in rechargeable backup battery allows continued use during power loss
- Includes two side pockets and a hand control hook for added user convenience
- Optional extras available, including headrest pillow and pressure care seating options

Specifications

- Seat Height: 50 cm
- Seat Depth: 52 cm
- Seat Width: 69 cm
- Arm Height: 19 cm
- Back Height: 72.5 cm
- Overall Height: 115 cm
- Overall Depth: 100 cm
- Safe Working Load: 370 kg