



# AMBASSADOR VINYL PREMIER A4 BARIATRIC RECLINER

\$6,976.00 - \$9,734.00Price range: \$6,976.00 through \$9,734.00



**SKU: PA4VINYL-H** 

#### **PRODUCT DESCRIPTION**

### **Purpose and Benefits**

The Ambassador Vinyl Premier A4 Bariatric Lift Chair is designed to provide maximum comfort, support, and independence for individuals who require a higher weight capacity seating solution. With a safe working load of up to 370 kg, this chair is ideal for users seeking a durable and reliable lift and recline chair that allows smooth and safe transitions between sitting and standing.

Its heavy-duty construction, supportive cushioning, and customisable design make it suitable for both home and healthcare environments. The A4 model combines functional assistive features with highend comfort, helping to improve daily living for users with mobility challenges.

#### **Features**

- Four-motor design allows independent operation of backrest and footrest
- Leg rest elevates 85 mm above seat height to improve lower limb circulation
- Reinforced steel frame construction for long-lasting structural support
- Customisable seat width from 69 cm to 95 cm to suit user needs
- Upholstery covers are removable and replaceable for hygiene and convenience
- Rechargeable backup battery enables operation during power outages



- Two built-in side pockets and a hand control hook for easy access
- Optional accessories include armrest covers, headrest pillow, and pressure care cushions

## **Specifications**

• Seat Height: 50 cm

• Seat Depth: 52 cm

• Seat Width: 69 cm to 95 cm (customisable)

• Arm Height: 19 cm

• Back Height: 72.5 cm

• Overall Height: 115 cm

• Overall Depth: 100 cm

• Safe Working Load: Up to 370 kg

• Upholstery Options: Soft-touch, high-quality fabric with removable covers

• Warranty: 5 years on frame, 3 years on electronics and motors