



AMBASSADOR VINYL PREMIER A4 BARIATRIC RECLINER



SKU: PA4VINYL-H

PRODUCT DESCRIPTION

Purpose and Benefits

The Ambassador Vinyl Premier A4 Bariatric Lift Chair is designed to provide maximum comfort, support, and independence for individuals who require a higher weight capacity seating solution. With a safe working load of up to 370 kg, this chair is ideal for users seeking a durable and reliable lift and recline chair that allows smooth and safe transitions between sitting and standing.

Its heavy-duty construction, supportive cushioning, and customisable design make it suitable for both home and healthcare environments. The A4 model combines functional assistive features with high-end comfort, helping to improve daily living for users with mobility challenges.

Features

- Four-motor design allows independent operation of backrest and footrest
- Leg rest elevates 85 mm above seat height to improve lower limb circulation
- Reinforced steel frame construction for long-lasting structural support
- Customisable seat width from 69 cm to 95 cm to suit user needs
- Upholstery covers are removable and replaceable for hygiene and convenience
- Rechargeable backup battery enables operation during power outages

- Two built-in side pockets and a hand control hook for easy access
- Optional accessories include armrest covers, headrest pillow, and pressure care cushions

Specifications

- Seat Height: 50 cm
- Seat Depth: 52 cm
- Seat Width: 69 cm to 95 cm (customisable)
- Arm Height: 19 cm
- Back Height: 72.5 cm
- Overall Height: 115 cm
- Overall Depth: 100 cm
- Safe Working Load: Up to 370 kg
- Upholstery Options: Soft-touch, high-quality fabric with removable covers
- Warranty: 5 years on frame, 3 years on electronics and motors

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au