




AQUASENSE MULTI ADJUST BATH SAFETY RAIL

\$87.00



SKU: 785-340

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The new AquaSense multi-adjust bath safety rail allows for both parallel and perpendicular orientation to the bathtub wall, and features four different height settings. It provides a feeling of security, stability and comfort when entering or exiting the bathtub.

FEATURES

- Heavy gauge, rust-resistant stainless steel and aluminum construction.
- Two independent contact surfaces with the tub wall, for strong and reliable support.
- Fits tub walls from 3" to 7" (7.6 cm to 17.8 cm) thick.
- Non-marring rubber pads on all points of contact with the tub, for a solid yet bathtub-friendly grip.
- Installs securely without the use of tools.
- Weight capacity: 300 lb (136 kg)