



AQUASENSE BATH SAFETY RAIL



SKU: 785-350

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Added support for getting in and out of the tub!

This AquaSense Safety Rail provides safety and comfort while getting in and out of the bathtub. Designed to fit a wide variety of tubs, it offers two gripping surfaces. It installs easily and securely without the use of tools and is manufactured from rust-resistant steel. It features an easy to clean, white vinyl coated surface, and non-marring rubber pads on all contact points.

FEATURES

- Non-marring rubber pads on all contact points for a solid yet bathtub-friendly grip.
- Large adjustment knob offers a good grip and allows for easy, tool-free installation.
- Two independent gripping surfaces allow this safety rail to adapt to a wide variety of tubs.
- Heavy-gauge, rust-resistant welded steel construction offers solid support.
- Easy to clean, white vinyl coated surface.
- Fits bathtub side walls from 3" to 7" (7.6 to 17.8 cm) wide.
- Height above bathtub: 15" (38 cm).
- Maximum capacity: 300 lb (136 kg)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au