



## **AQUATEC OCEAN ERGO SELF PROPELLED SHOWER COMMODE**

**~~\$2,933.00~~ Original price was:**

**~~\$2,933.00.~~ \$2,639.70 Current price is: \$2,639.70.**



**SKU:** PKG-1643718SP

## **PRODUCT DESCRIPTION**

### **PURPOSE AND BENEFITS**

The Aquatec Ocean Ergo Self Propelled is a high-quality shower chair designed for ease of use and maximum adjustability. Engineered by Aquatec in Germany, this chair boasts a stable stainless steel frame and tool-less height adjustable seat for a customized experience. The ergonomic swing-away armrests, height adjustable swing-away footrests, and fold-away texture footplates provide comfort and safety for both the user and carer. The machine washable backrest and space-saving push bar make this chair convenient and easy to use. With a maximum weight capacity of 150 kg (with armrest width extension), the Aquatec Ocean Ergo Self Propelled is a reliable and efficient solution for those who need assistance in the shower or toilet.

### **FEATURES**

- Stable Stainless Steel Frame
- Tool-less, height-adjustable seat and footplate
- Ergonomic, flip-up armrests are also height adjustable (Allen key required - supplied)
- Washable, tension-adjustable backrest with quick-release buckles
- Ergonomically-designed seat
- Four braked castors for optimum safety
- Easy to clean with removable seat and smooth surfaces
- Full range of accessories available (NB: Transit version can be retrofitted to self-propelled model with 610mm/24" wheels)
- Maximum weight capacity of 150Kg (with armrest width extension)
- Comes standard with Pan Holder and Soft Seat.

## SPECIFICATIONS

Seat Height:	47 to 60cms
Seat Width:	48cm
Seat Depth:	45cm
Overall Width:	67.5cm
Overall Length:	106cm
Length without Footrests:	73cm
Footrest Depth:	15.4cm
Overall Height:	96.4 to 109cms
Height to Armrests:	68 to 81cms
Height - Footrests to Seat:	28 to 53cms
Height: ground to lower seat edge with holder for pan:	39 to 52cms
Height: ground to lower seat edge without holder for pan:	41 to 54cms
Width between Armrests:	46cm (54cm with XL extension upgrade kit)
Height Adjustable Steps:	2.5cm
Overall Weight:	19.5kg
Maximum Weight Capacity:	150kg