



## AQUATEC OCEAN DUAL VIP ERGO SHOWER COMMODE



**SKU:** PKG-1643717

## PRODUCT DESCRIPTION

### Purpose & Benefits of the Aquatec Ocean Dual VIP Ergo Shower Commode

The Aquatec Ocean Dual VIP Ergo Shower Commode is expertly crafted to provide exceptional comfort and versatility for individuals with limited mobility. This advanced shower commode chair features both tilt-in-space and recline functionalities that enable precise positioning to accommodate various therapeutic needs, enhance comfort, and assist in pressure management. It is especially useful for those requiring support during hygiene routines and is ideal for individuals at risk of pressure sores, or who need assistance maintaining posture due to medical conditions.

### Features

- **Dual Tilt and Recline Mechanisms:** Offers independent tilt and recline options that allow for optimal positioning to reduce pressure risk, improve comfort, and support complex care needs.
- **Ergonomic Seating:** The seat design conforms to the natural shape of the user, providing anatomical support. Adjustable headrests and footrests further enhance comfort and accommodate varying body sizes.
- **Modular and Adjustable Design:** The chair's modular design allows for numerous adjustments including seat height, footrest position, and armrest width, making it highly adaptable to individual user needs.
- **Durable, Corrosion-Resistant Materials:** Built to withstand the humid conditions of shower environments, using corrosion-resistant materials to ensure durability and longevity.
- **Safety Features:** Includes lockable casters to secure the chair in place during transfers and use, as well as safety belts for additional security.
- **Ease of Cleaning:** Smooth surfaces and removable components simplify the cleaning process,

essential for maintaining hygiene and extending the product's lifespan.

## Specifications

- **Material:** Constructed from high-quality stainless steel and durable, water-resistant plastics.
- **Safe Working Load:** 150kg
- **Adjustability:** Features extensive adjustability in seat height, backrest angle, and tilt-in-space mechanism.
- **Dimensions:**
  - **Seat Height:** 51 - 61cm
  - **Seat Width:** 50cm
  - **Total Width:** 59.5cm
  - **Total Depth:** 106cm
  - **Total Height without Headrest:** 108.5 - 118.5cm
- **Care Instructions:** Regular washing with mild, non-abrasive cleaners is recommended. Check and maintain all moving parts periodically to ensure optimal function and safety.

The Aquatec Ocean Dual VIP Ergo Shower Commode addresses the complex needs of users with significant mobility limitations, offering a comprehensive solution that enhances their quality of life. Its ability to adjust to various care requirements and settings makes it an invaluable asset in both home care and medical facilities, promoting independence, comfort, and dignity for users.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### Monash Showroom

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)