



AQUATEC OCEAN ERGO TRANSIT COMMODE

\$1,812.00



SKU: PKG-1643718

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Introducing the Ocean Ergo Transit Commode, a thoughtfully designed aid for individuals with mobility challenges. This commode chair is engineered to provide a comfortable and safe bathroom experience, with features that cater to both the user's and caregiver's needs.

This kit also comes with a **Pan Holder** and **Soft Seat**.

FEATURES

- **Ergonomic Seat:** The ergonomically designed seat enhances comfort and accessibility, making personal hygiene tasks easier and more dignified.
- **Adjustable Footrests:** Footrests can be adjusted in height and angle to accommodate different user requirements, promoting proper posture and support.
- **Swing-Away Armrests:** For ease of transfer and to assist caregivers, the armrests swing away smoothly, providing unobstructed access.
- **Transit Wheels:** Equipped with four transit wheels, the commode can be easily maneuvered over toilets and into showers, with locks to ensure stability during use.
- **Corrosion-Resistant Frame:** Constructed with high-quality, corrosion-resistant materials suitable for use in wet environments.
- **Safety Features:** Includes brakes on all wheels and a stable base to ensure user safety while sitting or transferring.

The Ocean Ergo Transit Commode stands out for its user-centric design, providing a safe and convenient solution for those who need assistance with daily toileting and showering tasks. Its blend

of comfort, safety features, and ease of use makes it a valuable addition to any home or healthcare facility.

SPECIFICATIONS

Overall Width	565 mm
Overall Depth	905 mm
Seat Width	480 mm
Seat Depth	480 mm
Min Seat Height	475 mm
Max Seat Height	600 mm
Product Weight	17 kg
Max User Weight	150 kg
Brand	KCare