



## AQUATEC OCEAN ERGO TRANSIT COMMODE



**SKU:** PKG-1643718

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

Introducing the Ocean Ergo Transit Commode, a thoughtfully designed aid for individuals with mobility challenges. This commode chair is engineered to provide a comfortable and safe bathroom experience, with features that cater to both the user's and caregiver's needs.

This kit also comes with a **Pan Holder** and **Soft Seat**.

### FEATURES

- **Ergonomic Seat:** The ergonomically designed seat enhances comfort and accessibility, making personal hygiene tasks easier and more dignified.
- **Adjustable Footrests:** Footrests can be adjusted in height and angle to accommodate different user requirements, promoting proper posture and support.
- **Swing-Away Armrests:** For ease of transfer and to assist caregivers, the armrests swing away smoothly, providing unobstructed access.
- **Transit Wheels:** Equipped with four transit wheels, the commode can be easily maneuvered over toilets and into showers, with locks to ensure stability during use.
- **Corrosion-Resistant Frame:** Constructed with high-quality, corrosion-resistant materials suitable for use in wet environments.
- **Safety Features:** Includes brakes on all wheels and a stable base to ensure user safety while sitting or transferring.

The Ocean Ergo Transit Commode stands out for its user-centric design, providing a safe and convenient solution for those who need assistance with daily toileting and showering tasks. Its blend

of comfort, safety features, and ease of use makes it a valuable addition to any home or healthcare facility.

## SPECIFICATIONS

Overall Width	565 mm
Overall Depth	905 mm
Seat Width	480 mm
Seat Depth	480 mm
Min Seat Height	475 mm
Max Seat Height	600 mm
Product Weight	17 kg
Max User Weight	150 kg
Brand	KCare

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
 Brighton SA 5048  
 08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
 Ave, Elizabeth Vale SA 5112  
 08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,  
 Reynella SA 5161  
 08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199  
 03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### Monash Showroom

23 Winterton Road, Clayton VIC 3168  
 03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

**Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)