




AQUATEC OCEAN VIP ERGO SHOWER COMMODE



SKU: PKG-1643577

PRODUCT DESCRIPTION

Purpose & Benefits of Aquatec Ocean VIP Ergo

The Aquatec Ocean VIP Ergo is a versatile and highly adaptive shower commode chair designed to provide enhanced comfort and support for individuals with limited mobility. It features a tilt-in-space mechanism that helps manage pressure redistribution, which is crucial for users at risk of skin breakdown or those who require prolonged periods of sitting. The VIP Ergo model is especially beneficial for users needing additional postural support and those undergoing rehabilitation or dealing with severe mobility impairments.

Features

- **Tilt-in-Space Functionality:** Allows the chair to tilt up to 35 degrees, providing optimal pressure relief and comfort. This feature aids in the prevention of pressure sores and enhances physiological functions.
- **Ergonomic Seat Design:** Designed with contoured seating that conforms to the body's natural shape, reducing pressure points and enhancing user comfort during use.
- **Adjustable Components:** Includes adjustable headrests, footrests, and armrests, allowing for a fully customized fit that meets the specific needs of each user.
- **High-Quality Materials:** Constructed from durable, corrosion-resistant materials that are ideal for use in wet environments like bathrooms.
- **Safety Features:** Equipped with lockable casters to ensure stability and safety when in use. The chair also includes a safety belt for added security during transfers or while showering.
- **Easy to Clean:** Features smooth surfaces and removable cushions that make it easy to clean and maintain, ensuring high levels of hygiene are upheld.

Specifications

- **Material:** High-grade stainless steel and durable plastic components, which are corrosion-resistant and suitable for wet environments.
- **Weight Capacity:** Capable of supporting a variety of weights, typically accommodating up to 150 kg (330 lbs).
- **Adjustability:** Extensive adjustability features including tilt-in-space, backrest angle, and footrest height.
- **Dimensions:**
 - Width: 56.5cm
 - Length: 100cm
 - Product Weight: 23.5kg
 - Safe Working Load: 150kg
 - Seat Height: 47.5cm - 60cm
 - Seat Width: 48cm
 - Seat Depth: 48cm
 - Seat Tilt: -5° to 40°
- **Care Instructions:** Easy to maintain by regularly cleaning with mild soap and water, especially after use in wet conditions.

The Aquatec Ocean VIP Ergo shower commode chair is an invaluable tool for enhancing the quality of life for individuals with significant mobility challenges. Its ergonomic design and adjustable features ensure that users can maintain hygiene independently or with minimal assistance while receiving the necessary support to sit comfortably and safely during use.

For further information go to **[medimart.com.au](https://www.medimart.com.au)** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au