



AQUATEC ORCA BATH LIFT



SKU: 1470753-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Aquatec Orca Bath Lift is a safe and comfortable solution designed to help individuals with limited mobility to enjoy the pleasure of a bath. With its reclining backrest and sturdy design, it offers both safety and comfort.

FEATURES

- **Reclining Backrest:** Provides a comfortable bathing position while ensuring safety and stability.
- **Easy-to-Use Controls:** The waterproof, floating hand control is simple to use, making it easy for individuals of all ages to operate.
- **Battery Powered:** Typically equipped with a rechargeable battery, ensuring that the lift can be operated without the need for a continuous electrical supply, enhancing safety in a wet environment.
- **Safety Features:** Many models come with built-in safety mechanisms that prevent the lift from lowering if there isn't enough battery power to raise it back up, ensuring that users are not left stranded in the tub.
- **Comfortable Seating:** The high-quality materials used for the seating area provide comfort and are easy to clean.
- **Portable Design:** In some cases, bath lifts are designed to be lightweight and easy to disassemble, making them portable and easy to store.
- **Non-Marking Suction Feet:** Secure the lift firmly to the bathtub, ensuring it stays in place during use, and do not leave marks on the tub.

The Aquatec Orca Bath Lift is an innovative product that enhances the bathing experience for

individuals facing challenges with mobility. It combines functionality with ease of use, offering a secure and relaxing bath time.

SPECIFICATIONS

Aquatec Orca Standard

Overall Width	715 mm
Overall Depth	650 mm
Overall Height	1075 mm
Min Seat Width	375 mm
Max Seat Width	715 mm
Seat Depth	590 mm
Min Seat Height	60 mm
Max Seat Height	420 mm
Backrest Height	665 mm
Product Weight	10.8 kg
Max User Weight	140 kg

Aquatec Orca Extra Large

Overall Width	715 mm
Overall Depth	650 mm
Overall Height	1075 mm
Min Seat Width	375 mm
Max Seat Width	715 mm
Seat Depth	590 mm
Min Seat Height	60 mm
Max Seat Height	420 mm
Backrest Height	665 mm
Product Weight	11.6 kg
Max User Weight	170 kg

<https://youtu.be/dBg9pAlGoSo>

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au