


NOVIS AREACARE MATTRESS EXTENSION

\$123.00 - \$147.00



SKU: BSPS5400-H



PRODUCT DESCRIPTION

Purpose and Benefits

This clinically designed mattress extension adds 200mm of extra length to compatible pressure care mattresses, offering enhanced support for taller individuals or those requiring focused heel protection. It helps maintain proper body alignment while delivering effective pressure redistribution in critical areas such as the heels, reducing the risk of pressure injuries.

The dual-layer foam construction features a high-density memory foam top for targeted comfort and a medium-density core for stability and durability. A waterproof, multi-stretch cover with welded and taped seams enhances hygiene and infection control. Ideal for use in hospitals, aged care facilities, and home environments, the extension integrates seamlessly with selected 150mm-high mattresses—extending comfort and functionality without replacing the full mattress.

Specifications

Single

- Length: 200 mm
- Width: 865 mm
- Height: 150 mm
- Weight capacity: As per accompanying mattress used

King Single

- Length: 200 mm
- Width: 1050 mm
- Height: 180 mm
- Weight capacity: As per accompanying mattress used

Warranty

- Cover: 2 years
- Body: 6 years

Compliance

- ARTG Number: 176837
- Australian Standards: BS7177, BS5852