


NOVIS AREACARE MATTRESS EXTENSION



SKU: BSPS5400-H

PRODUCT DESCRIPTION

Purpose and Benefits

This clinically designed mattress extension adds 200mm of extra length to compatible pressure care mattresses, offering enhanced support for taller individuals or those requiring focused heel protection. It helps maintain proper body alignment while delivering effective pressure redistribution in critical areas such as the heels, reducing the risk of pressure injuries.

The dual-layer foam construction features a high-density memory foam top for targeted comfort and a medium-density core for stability and durability. A waterproof, multi-stretch cover with welded and taped seams enhances hygiene and infection control. Ideal for use in hospitals, aged care facilities, and home environments, the extension integrates seamlessly with selected 150mm-high mattresses—extending comfort and functionality without replacing the full mattress.

Specifications

Single

- Length: 200 mm
- Width: 865 mm
- Height: 150 mm
- Weight capacity: As per accompanying mattress used

King Single

- Length: 200 mm
- Width: 1050 mm
- Height: 180 mm
- Weight capacity: As per accompanying mattress used

Warranty

- Cover: 2 years
- Body: 6 years

Compliance

- ARTG Number: 176837
- Australian Standards: BS7177, BS5852

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA**Brighton Showroom**

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC**Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au