



ARJO 4 POINT DOUBLE AMPUTEE SLING WITH HEAD SUPPORT



SKU: MAA2090M-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Suitable for patients with bilateral leg amputations above the knee or total hip, this sling provides good comfort and support for a range of everyday patient handling routines, such as bed and chair-related transfers.

Designed for use with four-point Dynamic Positioning System (DPS) spreader bars, which provide a more comfortable reclined position in the sling and a more even distribution of the resident's/patient's weight.

The DPS system allows the caregiver to position the resident/patient in the sling from upright to reclined with minimum physical effort.

FEATURES

- Suitable for patients with bilateral leg amputations above the knee or total hip.
- Designed to ensure good comfort and support for a range of everyday patient handling routines, such as bed and chair-related transfers.
- Made up of a soft and durable polyester fabric.
- DPS system allows the caregiver to position the patient in the sling.
- Available in different sizes: XS, S, M, L, XL, XXL

SPECIFICATIONS

- Safe working load: 272 kg/600 lbs

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au