



## **ARJO 4 POINT DOUBLE AMPUTEE SLING WITH HEAD SUPPORT**

**\$647.00**



**SKU:** MAA2090M-H

## **PRODUCT DESCRIPTION**

### **PURPOSE & BENEFITS**

Suitable for patients with bilateral leg amputations above the knee or total hip, this sling provides good comfort and support for a range of everyday patient handling routines, such as bed and chair-related transfers.

Designed for use with four-point Dynamic Positioning System (DPS) spreader bars, which provide a more comfortable reclined position in the sling and a more even distribution of the resident's/patient's weight.

The DPS system allows the caregiver to position the resident/patient in the sling from upright to reclined with minimum physical effort.

### **FEATURES**

- Suitable for patients with bilateral leg amputations above the knee or total hip.
- Designed to ensure good comfort and support for a range of everyday patient handling routines, such as bed and chair-related transfers.
- Made up of a soft and durable polyester fabric.
- DPS system allows the caregiver to position the patient in the sling.
- Available in different sizes: XS, S, M, L, XL, XXL

### **SPECIFICATIONS**

- Safe working load: 272 kg/600 lbs