



ARJO 4 POINT MESH BATHING SLING WITH HEAD SUPPORT



SKU: MAA4060M-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Mesh material allows water to pass easily through the sling, making it an ideal choice for showering or bathing. The DPS system allows the caregiver to position the resident/patient in the sling from upright to reclined with minimum physical effort. With no padding in the leg area, the unpadded sling is easier to apply to residents/patients with skeletal or spinal conditions.

FEATURES

- Made up of a soft and durable polyester fabric.
- DPS system allows the caregiver to position the resident/patient in the sling.
- Suitable for handling routines, such as bed and chair-related transfers.
- Maximum patient comfort
- Minimum physical efforts of caretaker
- Available in different sizes: M, LL, XL

SPECIFICATIONS

- Safe working load: 272 kg/600 lbs

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au