



## **ARJO 4 POINT MESH BATHING SLING WITH HEAD SUPPORT**

**\$547.00 - \$563.00** Price range: \$547.00 through \$563.00



**SKU:** MAA4060M-H

## **PRODUCT DESCRIPTION**

### **PURPOSE & BENEFITS**

Mesh material allows water to pass easily through the sling, making it an ideal choice for showering or bathing. The DPS system allows the caregiver to position the resident/patient in the sling from upright to reclined with minimum physical effort. With no padding in the leg area, the unpadded sling is easier to apply to residents/patients with skeletal or spinal conditions.

### **FEATURES**

- Made up of a soft and durable polyester fabric.
- DPS system allows the caregiver to position the resident/patient in the sling.
- Suitable for handling routines, such as bed and chair-related transfers.
- Maximum patient comfort
- Minimum physical efforts of caretaker
- Available in different sizes: M, LL, XL

### **SPECIFICATIONS**

- Safe working load: 272 kg/600 lbs