



## **ARJO 4 POINT MESH BATHING SLING WITH HEAD SUPPORT**



**SKU:** MAA4060M-H

## **PRODUCT DESCRIPTION**

### **PURPOSE & BENEFITS**

Mesh material allows water to pass easily through the sling, making it an ideal choice for showering or bathing. The DPS system allows the caregiver to position the resident/patient in the sling from upright to reclined with minimum physical effort. With no padding in the leg area, the unpadded sling is easier to apply to residents/patients with skeletal or spinal conditions.

### **FEATURES**

- Made up of a soft and durable polyester fabric.
- DPS system allows the caregiver to position the resident/patient in the sling.
- Suitable for handling routines, such as bed and chair-related transfers.
- Maximum patient comfort
- Minimum physical efforts of caretaker
- Available in different sizes: M, LL, XL

### **SPECIFICATIONS**

- Safe working load: 272 kg/600 lbs

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)