



## **ARJO ALPHA ACTIVE 4 ALTERNATING PRESSURE SYSTEM**

**\$800.00 - \$1,387.00** Price range: \$800.00 through \$1,387.00



**SKU:** 648326-H

### **PRODUCT DESCRIPTION**

#### **PURPOSE AND BENEFITS**

Designed for the prevention and management of pressure injuries, the Alpha Active 4 pressure redistributing mattress replacement system is suitable for a range of care environments, including home and community care settings.

The Alpha Active 4 system has been designed as a cost-effective solution to aid in the prevention and management of pressure injuries.

The system offers advanced features for the delivery of high quality care to vulnerable patients and residents, whilst remaining easy-to-use for the carer and simple to clean and maintain.

The Arjo Alpha Active 4 system requires both the mattress replacement system and pump to be purchased.

#### **FEATURES**

- Electronic pressure setting allows accurate and appropriate pressure settings to be set for various weight ranges
- Variable cycle times
- AutoFirm provides a stable support surface for patient transfer and nursing procedures
- Two modes of operation: active / alternating and reactive / static
- Includes low pressure alarm, power fail alarm, and pump fault alarm
- Modular construction contributes to ease of cleaning and maintenance

- Ventilated standard air cells designed to aid airflow within the mattress
- Cell-on-Cell designed to provide additional support and to help increase patient stability
- Three static head cells
- Reliant IS2 cover is vapour permeable and water resistant, incorporating 2-way stretch and bacteriostatic properties
- Rapid deflation for CPR
- **Pump sold separately**
- **Both the Mattress Replacement System and Pump is required to operate the Arjo Alpha 4 Alternating Pressure System**

## SPECIFICATIONS

Total Length	Total Width	Total Height	SWL
1950mm	850mm	200mm	200kg