



ARJO ATMOS AIR VELARIS ALTERNATING PRESSURE MATTRESS

\$1,550.00 - \$3,299.00 Price range: \$1,550.00 through \$3,299.00



SKU: 633021EN-H

PRODUCT DESCRIPTION

Purpose & Benefits of the Arjo AtmosAir Velaris Alternating Pressure Mattress

The Arjo AtmosAir Velaris Alternating Pressure Mattress is an advanced pressure care solution designed for individuals at high risk of developing pressure injuries. Combining innovative self-adjusting air-cell technology with alternating pressure functionality, this mattress ensures optimal pressure redistribution and skin integrity. Ideal for use in clinical, aged care, and homecare environments, the AtmosAir Velaris promotes patient comfort, safety, and effective pressure injury prevention.

Features

- **Self-Adjusting Air Cells:** The mattress features self-adjusting air cells that automatically adapt to the user's weight and position, delivering consistent pressure redistribution and comfort.
- **Alternating Pressure Therapy:** Provides alternating pressure cycles to reduce pressure on key areas, enhancing circulation and reducing the risk of pressure injuries.
- **Quiet and Efficient Pump System:** Equipped with a low-noise pump for minimal disruption, the system ensures continuous operation and reliability.
- **Integrated Foam Layer:** Combines air-cell technology with a supportive foam layer for enhanced comfort and postural support.
- **Moisture-Resistant Cover:** The vapour-permeable, waterproof cover protects the mattress from spills and promotes hygiene while supporting infection control.
- **Customisable Settings:** Adjustable settings allow caregivers to tailor the therapy to meet individual patient needs.
- **Durable and Long-Lasting:** Constructed with high-quality materials to ensure reliable performance in demanding care settings.

- **THE PUMP IS SOLD SEPARATELY AND MUST BE PURCHASED IN ADDITION TO THE MATTRESS.**

Specifications

AtmosAir Velaris Range

Top cover fabric options	Dartex Reliant or Premium
PVC free	Yes
Welded seams & handles as standard	Yes
Mattress fire tests***	BS7177 (BS597-1, BS597-2, BS6807 Clause 9-Ignition Source 5), 16CFR1632, 16CFR1633, CAL TB 129
Firesock as standard	Yes

AtmosAir Velaris Standard Mattress

Single	860 x 1980 x 180 mm (34 x 78 x 7 in)
King Single	880 x 2020 x 180 mm (35 x 80 x 7 in)
Ability to connect to AtmosAir Velaris alternating pressure pump	Yes
Weight of mattress	15.0 kg (33 lb) – 17.5 kg (39 lb)
Maximum patient weight	250 kg (550 lb)
Expected service life	5 years

AtmosAir Velaris Alternating Pressure Pump

Size (W x L x H)	337 x 107 x 200 mm (13.3 x 4.2 x 7.9 in)
Weight	4.1 kg (9 lb)
Alternating cycle time	10 minutes
Supply voltage	100 – 230 V
Power input	3 – 46 VA
IP rating	IP22
Case material	PC ABS
Sound output	19 dBA
Lock out control panel feature	Yes
Suppress audible Alarm option	Yes
Alarm pause	Yes (15 mins)
Skin IQ adapter part number**	Power cable 636377
Expected service life	7 years

The Arjo AtmosAir Velaris Alternating Pressure Mattress is a cutting-edge solution for pressure injury prevention, combining comfort, durability, and innovative technology to deliver superior care and support.