



ARJO AURALIS STANDARD MATTRESS SINGLE 200



SKU: 636M02S

PRODUCT DESCRIPTION

Purpose & Benefits

The **Arjo Auralis Standard Mattress Single 200** is an air-filled zone sub-mattress designed exclusively for use with the Arjo Auralis Alternating Pressure Mattress System. It is intended for high acuity patients with limited mobility and compromised skin integrity, helping to minimise the risk and burden of pressure injuries. The mattress supports both alternating pressure and constant low pressure therapy through the Auralis pump's Self-Set Technology, which automatically adjusts pressure levels to suit individual patient needs. Additional features including Micro Air-Loss, Heelguard® technology, Cell in Cell construction, transport mode, CPR control, and a breathable water-resistant cover are designed to enhance patient comfort, stability, hygiene, and clinical care.

Features

- Compatible with the Arjo Auralis Alternating Pressure Mattress System
- Designed for patients with limited mobility and compromised skin integrity
- Supports alternating pressure and constant low pressure therapy
- Self-Set Technology automatically adjusts mattress pressure
- Micro Air-Loss function helps reduce heat and moisture build-up
- Heelguard® technology assists with heel pressure redistribution
- Cell in Cell technology improves stability during transfers

- Transport mode supports the patient for up to 12 hours
- Water-resistant breathable cover with welded seams
- CPR control enables rapid mattress deflation
- Integrated cable management reduces trailing cables

Specifications

Model	636M02S
Mattress Size	203 cm × 86 cm × 20.5 cm
Mattress Type	Air-filled zone sub-mattress
Compatibility	Arjo Auralis Alternating Pressure Mattress System
Cover	Water-resistant breathable premium cover with welded seams

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048

08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112

08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriiffs Rd,
Reynella SA 5161

08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au