



ARJO BARIATRIC LOOP SLING

\$1,070.00



SKU: VIG11004-L1

PRODUCT DESCRIPTION

Purpose & Benefits

The **Arjo Bariatric Loop Sling** is designed to assist with transfers of plus-size patients or residents who require full-body support and have limited ability to move. It enables safe transfers using Arjo bariatric passive lifts, reducing risk of injury for both users and caregivers, and improving comfort through padded leg pieces.

Features

- Standard padding in the legs to enhance comfort
- Compatible with Arjo passive lift systems using loop-attachment only
- Swings between bed, chair, wheelchair, bathroom and floor in seated position
- Meets EN ISO 10535 standard for hoist slings

Specifications

- **Product Code:** MAA8010
- **Safe Working Load:** 454 kg

- **Compliance:** EN ISO 10535
- **Warranty:** 3 months (manufacturer's defects)