



ARJO CALYPSO HEIGHT-ADJUSTABLE BATH AND SHOWER CHAIR



SKU: CDB6004-01

PRODUCT DESCRIPTION

Purpose & Benefits of the Arjo Calypso Height-Adjustable Bath and Shower Chair

The Arjo Calypso is a versatile, height-adjustable bath and shower chair designed to facilitate safe and efficient hygiene routines for both caregivers and patients. It enables a single caregiver to transfer patients directly from their bedside, chair, or wheelchair to the bath, shower, or toilet without manual lifting, thereby reducing the risk of injury and enhancing comfort. The Calypso is compatible with all height-adjustable baths from Arjo, making it a flexible solution for various care settings.

Features

- **Height Adjustability:** The chair can be raised or lowered to align with different surfaces, ensuring smooth transfers and accommodating various bath heights.
- **Foldable Backrest and Armrests:** Individually foldable back and armrests allow for easy patient positioning and comfort during bathing or showering.
- **Ergonomic Design:** Equipped with easy-to-grip handles for manoeuvring and positioning, enhancing caregiver control and patient security.
- **Versatile Seat Options:** Available with a classic seat or a soft, easy-to-clean polyurethane material, both featuring hygiene openings to facilitate toileting.
- **Safety Features:** Includes a manually applied foot brake on all castors to ensure stability during transfers and bathing routines.
- **Optional Accessories:** Additional features such as a bedpan with holder and a foot scale are available to meet specific patient needs.

Specifications

- **Maximum Safe Working Load:** 136 kg
- **Lifting Stroke:** 500 mm
- **Total Width:** 631 mm
- **Seat Width:** 479 mm
- **Total Length:** 933 mm
- **Seat Height Range:** 470 mm to 970 mm
- **Unit Weight:** 55 kg

The Arjo Calypso Height-Adjustable Bath and Shower Chair is an essential tool for enhancing patient hygiene care, promoting dignity and comfort while ensuring caregiver safety and efficiency.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au