



ARJO CARENDO ELECTRIC HYGIENE SYSTEM

\$18,310.00



SKU: BIB2004-01

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The Carendo multipurpose hygiene chair facilitates an efficient hygiene and showering routine within a single transfer, while allowing a comfortable working posture for the caregiver throughout the process.

Carendo is designed to enable a single caregiver to perform the full hygiene routine, including dressing and undressing, toileting, showering and other hygiene tasks in a single transfer.

Height-adjustable and reclining, Carendo allows the caregiver and resident to reach their preferred height and position throughout the hygiene routine.

FEATURES

- Care Raiser enables the carer to gently lift the user's lower body with the touch of a button allowing easy access for performing discreet hygiene tasks
- Height adjustable and reclining
- Full length backrest
- Handheld control unit
- Emergency stop
- Emergency lowering
- Sturdy armrests
- Purchase includes: seat cushion | back cushion | two batteries | charging unit

SPECIFICATIONS

Minimum Length Lowest Position Upright	875mm
Maximum Length Highest Position Reclined	1510mm
Total Width	640mm
Chassis Width	640mm
Chassis Width between Legs	508mm
Chassis Length	865mm
Minimum Height Backrest Lowest Position Reclined	1020mm
Maximum Height Backrest Top Position Upright	1550mm
Minimum Height Armrests Lowest Position Reclined	775mm
Maximum Height Armrests Top Position Upright	1320mm
Seat Height Lowest Position	510mm
Seat Height Highest Position	1050mm
Lifting Stroke Seat	550mm
Weight	64kg
SWL	136kg