

ARJO CAREVO - HEIGHT ADJUSTABLE SHOWER TROLLEY

\$12,659.00



SKU: BAC1113-01



PRODUCT DESCRIPTION

ARJO CAREVO PURPOSE & BENEFITS

The Arjo Carevo shower trolley is designed to provide a comfortable shower experience for patients and residents who need to remain in a lying position, while enabling the caregiver to work efficiently and ergonomically.

Carevo is height-adjustable and features patented Ergo-Access areas to enable the caregiver to get closer to the resident and maintain a comfortable working posture, while the patented Flexi Zone is designed to enhance comfort for the patient or resident.

FEATURES

- Height adjustable shower trolley
- Facilitates improved working posture
- Clearly marked control panels located on both sides and at the head of the trolley, allowing easy access for the carer
- The soft surface material and patented Flexi Zone are designed to enhance patient or resident comfort throughout the hygiene routine
- Comfort handles provide support for the resident or patient
- Purchase includes: head pillow | battery (x2) | battery charger | battery holder | drainage hose – upper drainage | lower drainage
- Hire includes: head pillow | one battery | battery charger | battery holder | drainage hose – upper drainage | lower drainage

SPECIFICATIONS

Total Length	2005mm
Width Stretcher	586mm
Total Width Including Side Supports	759mm - 888mm
Chassis Length	1211mm
Chassis Width	700mm
Height Top Position Floor to Stretcher	955mm
Height Low Position Floor to Stretcher	605mm
Stroke Range	350mm
Weight Including Mattress	82kg
Maximum Patient Weight	182kg
SWL Patient and Water	200kg