


ARJO CAREVO - HEIGHT ADJUSTABLE SHOWER TROLLEY



SKU: BAC1113-01

PRODUCT DESCRIPTION

ARJO CAREVO PURPOSE & BENEFITS

The Arjo Carevo shower trolley is designed to provide a comfortable shower experience for patients and residents who need to remain in a lying position, while enabling the caregiver to work efficiently and ergonomically.

Carevo is height-adjustable and features patented Ergo-Access areas to enable the caregiver to get closer to the resident and maintain a comfortable working posture, while the patented Flexi Zone is designed to enhance comfort for the patient or resident.

FEATURES

- Height adjustable shower trolley
- Facilitates improved working posture
- Clearly marked control panels located on both sides and at the head of the trolley, allowing easy access for the carer
- The soft surface material and patented Flexi Zone are designed to enhance patient or resident comfort throughout the hygiene routine
- Comfort handles provide support for the resident or patient
- Purchase includes: head pillow | battery (x2) | battery charger | battery holder | drainage hose – upper drainage | lower drainage
- Hire includes: head pillow | one battery | battery charger | battery holder | drainage hose – upper drainage | lower drainage

SPECIFICATIONS

Total Length	2005mm
Width Stretcher	586mm
Total Width Including Side Supports	759mm - 888mm
Chassis Length	1211mm
Chassis Width	700mm
Height Top Position Floor to Stretcher	955mm
Height Low Position Floor to Stretcher	605mm
Stroke Range	350mm
Weight Including Mattress	82kg
Maximum Patient Weight	182kg
SWL Patient and Water	200kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
 Brighton SA 5048
 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
 Ave, Elizabeth Vale SA 5112
 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
 Reynella SA 5161
 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
 03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
 03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
 1300 633 468

support@medimarnt.com.au

