




## **ARJO CARINO - ADVANCED SHOWER SYSTEM**

**\$9,374.00**



**SKU: BOC1004-01**

## **PRODUCT DESCRIPTION**

### **ARJO CARINO PURPOSE AND BENEFITS**

Height-adjustable shower and hygiene chair

Designed to reduce the number of transfers required for assisted showering and hygiene, Carino enables a single caregiver to perform everyday hygiene tasks such as washing, toileting and foot care.

With Carino, showering and other hygiene routines can be carried out by one caregiver, thus facilitating personal one-to-one interactions.

## **FEATURES**

- Height-adjustable - Height-adjustability combined with a patented curved frame provide comfortable positioning and support for the patient or resident, while allowing the carer to reach their preferred working height at all times.
- Facilitates efficient showering routines - The patient or resident may remain seated in the Carino as they are transferred to and from the bedroom as well as for showering, toileting and other hygiene tasks, thus reducing the number of transfers required for daily hygiene routines.
- Swing-away leg rests and rotating footplates - Patented leg rests provide adjustable support when resident's legs are raised for hygiene tasks on the legs and feet, and can swing away to facilitate transfers, allow the resident to use the foot rests, or provide caregiver access to the patient's lower body. When not in use, the foot rests can also be tucked away under the Carino to allow the caregiver to move freely around the device.
- Suitable for mobility levels B and C.

## SPECIFICATIONS

Width mm	624 mm
Width in	24,7 in
Weight kg	42 kg
Weight lb	93 lb
Total width (incl. side supports) mm	630 mm
Total width (incl. side supports) in	24,75 in
Lifting stroke (hydraulic) mm	600 mm
Lifting stroke (hydraulic) in	23,7 in
Max Height mm	1310 mm
Max Height in	51,7 in
Low Height mm	881 mm
Low Height in	34,7 in
Safe Working Load kg	140 kg
Safe Working Load lb	309 lb
Total length of chassis mm	857 mm
Total length of chassis in	33,75 in
Width of chassis (between legs) mm	424 mm
Width of chassis (between legs) in	16,7 in
Backrest Max Height mm	1555 mm
Backrest Max Height in	61,25 in
Backrest Min Height mm	1240 mm
Backrest Min Height in	48,9 in
Seat Height low position mm	450 mm
Seat Height low position in	17,75 in
Seat Height High position mm	1050 mm
Seat Height High position in	41,4 in