




ARJO CARMINA - BARIATRIC SHOWER COMMODE



SKU: BLB2000-013

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

Designed for assisted showering, toileting and hygiene care for plus size patients and residents, Carmina offers a solution for transferring from the bedside, chair or wheelchair to the shower or toilet, and back again.

Carmina is suitable for plus size patients and residents with the mobility level of Albert, Barbara or Carl.

A solution for regular showering

Carmina is a mobile shower chair that supports the hygiene workflow, providing secure transport to and from the shower room and enabling hygiene care in a seated position.

Toileting

Carmina is designed to support the patient or resident during toileting routines by providing a commode function when access to the toilet is not possible.

Comfort for patient or resident and caregiver

Designed to accommodate the various body shapes and weight distributions of plus size patients, with comfort and dignity in mind. The adjustable backrest with self-adjusting angle allows a comfortable position to be achieved. Foldable armrests allow good access for the caregiver, whilst the additional push poles support maneuverability.

FEATURES

- Armrest - The armrests fold up, giving access for caregiver and enables lateral transfers.
- Adjustable backrest - Easily adjustable backrest has four positions and self-adjusting angle.
- Push poles - The push poles assist the caregiver when manoeuvring the chair.
- Tip-up seat - The tip-up seat is designed to make cleaning and disinfection easier and provides access to the commode container.

- Self-adjusting commode container holder - Container holder lifts up, self-adjusting to different toilet heights.
- Adjustable footrest - Pedal-operated footrest.

SPECIFICATIONS

Maximum Patient Weight, Safe Working Load kg	320 kg
Maximum Patient Weight, Safe Working Load lbs	705 lb
Total weight, excl. packaging (Carmina Basic) kg	37 kg
Total weight, excl. packaging (Carmina Basic) lbs	81 lb
Total weight, excl. packaging (Carmina) kg	75 kg
Total weight, excl. packaging (Carmina) lbs	165 lb
Maximum total weight (lift + patient) Carmina Basic kg	357 kg
Maximum total weight (lift + patient) Carmina Basic lbs	785 lb
Maximum total weight (lift + patient) Carmina kg	395 kg
Maximum total weight (lift + patient) Carmina lbs	869 lb
Overall width of product (Carmina Basic) mm	960 mm
Overall width of product (Carmina Basic) in	37 3/4 in
Overall width of product (Carmina) mm	997 mm
Overall width of product (Carmina) in	39 1/4 in
Overall depth of product (Carmina Basic) mm	760 mm
Overall depth of product (Carmina Basic) in	29 7/8 in
Overall depth of product (Carmina) mm	980 mm
Overall depth of product (Carmina) in	38 5/8 in
Overall height of product kg	1265 mm
Overall height of product lbs	49 3/4 in
Seat depths (4 positions) mm	420 to 525 mm
Seat depths (4 positions) in	16 1/2 to 20 5/8 in
Seat width mm	915 mm
Seat width in	36 in
Turning diameter mm	1300 mm
Turning diameter in	51 1/8 in
Ambient temperature range	+10 C to +40 C (+50 F to +104 F)
Water temperature range	+10 C to +40 C (+50 F to +104 F)
Relative humidity range	30% to 75%
Operating force brakes	<300 N

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au