

## ARJO CITADEL PLUS BARIATRIC BED



**SKU:** FXX21A4A4AKFBB

## PRODUCT DESCRIPTION

### Purpose & Benefits

The **Arjo Citadel Plus Bariatric Bed** is an advanced care system designed to provide safe, comfortable and efficient support for bariatric patients. With a Safe Working Load of up to 454 kg, this bed integrates intelligent features such as in-bed weighing, adjustable width, and safety monitoring to enhance both patient care and caregiver efficiency. It is ideal for hospital and acute care environments where precision, safety and dignity are priorities.

### Features

- **Integrated Weigh Scale:** Digital display accurately weighs up to 454 kg with  $\pm 500$  g precision. The patient's weight is displayed for 10 seconds to preserve dignity.
- **Variable Width & Length:** Bed frame expands to three adjustable width options to suit different patient sizes and care needs.
- **SafeSet™ Alerts:** Provides visual alerts for side rail position, brake status, bed height and VariZone activity, enhancing safety and awareness.
- **Compatible Support Surfaces:** Works with Arjo pressure redistribution systems for enhanced skin protection and comfort.
- **Head of Bed Angle Display:** Digital readout with automatic pause at 30° to support respiratory care and comfort.

- **VariZone™ Patient Movement Detection:** Detects and alarms when undesired patient movement or bed exit occurs.
- **AES (Anti-Entrapment System):** Automatically detects obstructions beneath the bed during articulation and halts downward movement.
- **One-Button Chair Position:** Effortlessly articulates the bed into a chair position to improve comfort and assist with patient mobility.

## Specifications

Specification	Details
Model	Arjo Citadel Plus Bariatric Bed
Safe Working Load (SWL)	454 kg
Weighing Accuracy	±500 g (digital integrated scale)
Width Adjustment	Three width expansion options
Head of Bed Angle	Digital display with 30° automatic pause
Safety System	SafeSet™ visual alerts; AES (Anti-Entrapment System); VariZone™ movement detection
Positioning	One-button chair position for enhanced comfort
Support Surface Compatibility	Compatible with Arjo pressure redistribution mattresses
Display Features	Digital weight and head angle display

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

**Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

**Monash Showroom**

23 Winterton Road, Clayton VIC 3168

03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

**NT**

**Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)