




## ARJO COMFORTFLX MATTRESS - 3 LAYERED FOAM



**SKU:** 1001246

### PRODUCT DESCRIPTION

### SIZE OPTIONS

Available in 8 different sizes - view the 'size table' tab above for specifications.

### PURPOSE AND BENEFITS

ComfortFLx foam mattresses have been designed and developed to provide patients and residents with a durable, high quality mattress, designed to meet the patients needs. - focusing on pressure injury prevention while offering a mattress comfortable enough for long term bed rest. The ComfortFLx is competitively priced, high quality multi layered hospital grade mattress, developed to immerse the patient in the mattress and provide effective pressure redistribution.

### KEY FEATURES

#### Cover Features

- Ultra two way stretch
- Vapour Permeable
- Water resistant
- Antimicrobial, bacteriostatic and fungistatic
- Concealed zip with welded seams
- Easy to clean and maintain
- Long lasting

#### Mattress features

- Patient weight of 0-250Kg
- Available in multiple sizes
- Mattress to fit your bedframe
- Ease of use
- Fire retardency to Australian standards

- High density hospital grade foam formula

#### Features and Benefits

##### Castellation

Castellations in the foam surface allow pressure redistribution and airflow for microclimate management, creating internal air channels to help reduce heat and moisture build-up while redistributing patient pressure for improved patient health and comfort.

##### Tri layer foam with firm side formers

Enables a safe working load of 250kg and a high quality pressure redistribution, meeting the requirements of the Australian hospital grade mattress.

##### Superb Profiling

Superb profiling in the gatched or supine position

Resident/patient protection and comfort The ComfortFLx range provides the optimum balance of effective pressure redistribution and a superior level of resident/patient comfort.

##### Unique design

Maximising the body contact area, the ComfortFLx range reduces the risks of localised areas of high pressure.