



ARJO ENTERPRISE 5000X - HOSPITAL BED



SKU: N/A

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The Enterprise 5000X hospital bed delivers significant benefits to the patient, caregiver and healthcare facility.

The daily challenges faced by nursing teams on hospital wards have inspired Arjo to improve upon an established solution, trusted by hospitals around the globe. In the development of the Enterprise 5000X hospital bed, Arjo has focused on key customer needs common to all healthcare environments, namely:

- Comfort for both resident/patient and caregiver.
- Compliance with new international standards.
- Equipment that is robust, intuitive and easy to use.
- Solutions focused on quality of care for residents/patients.

The bed can be supplied with or without 3/4 folding side rails.

FEATURES & BENEFITS

Security and comfort

A secure environment is critical to the wellbeing of patients and carers. Falls affect as many as 10% of older adults during an acute care inpatient stay¹ and can result in serious and often life-threatening injuries.

Intuitive operation and easy of use

Our one-button-push controls are intuitive and easy to use, thanks to clear and simple icons. One touch functions can be used for key features, such as:

- Auto CPR
- 30 degree backrest angle

- Auto-contour
- Cardiac chair position
- Low height and extra low height
- Trendelenberg / reverse Trendelenberg

Some or all of these functions may be locked out by the caregiver if required.

Enhancing the quality of care

The backrest will pause when raised to an angle of 30°, ensuring easy and positive backrest adjustment time after time. The caregiver can easily adjust the bed into a cardiac chair position with our simple to use control panel, by selecting the appropriate icon.

Infection control, pit stop maintenance and service

The Enterprise 5000X is designed with open architecture and no base cover to allow easy access for cleaning, decontamination, maintenance and service.

- The blow-moulded decksheets and head/foot panels are lightweight, smooth and easy to remove to aid cleaning.
- Pit stop maintenance is achieved by designing open accessibility and plug and play components. Changes such as actuator, control box and battery can be carried out in situ within a few minutes by one trained person.
- Equipment is designed to be easy to service and maintain, reducing time spent on cleaning, service and maintenance activities.

Optional: Now available with IndiGo

IndiGo is the powered 5th wheel that makes lighter transport available to any caregiver with just a touch – for ease and simplicity. IndiGo is always available and built directly into your Enterprise bed, so only one thing changes in a caregiver's day: the amount of work required to move the bed.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048

08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112

08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriiffs Rd,
Reynella SA 5161

08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au