


## ARJO EVALO - FLOORLINE BED



**SKU:** N/A

## PRODUCT DESCRIPTION

### PURPOSE AND BENEFITS

The EvaLo floor bed is designed to reduce the risk of harm from falls whilst providing you with a home-style bed for that creates peace of mind for you, your family and caregivers. Provides all the functions required for a comfortable sleeping experience and giving the confidence to move in the ways that you need. The ergonomic design allows it to fit in with practically any space with its clean design.

### FEATURES

#### Vertical High-Low

- Very low height – 160mm
- Electric Backrest, Knee bend and Height adjustment
- Trendelenburg and reverse Trendelenburg
- Integrated base extension
- Standard Head and Foot Boards – choice of colours
- SWL 200kg
- Available in King single