



## **ARJO GENERAL PURPOSE POLYESTER SLING WITH HEAD SUPPORT**

**\$500.00**



**SKU: MLA2000-H**

## **PRODUCT DESCRIPTION**

### **PURPOSE & BENEFITS**

Arjo standard padded loop sling is used on 2-point loop spreader bars. The slings are intended for everyday patient handling routines. The padded sling with head support offers padding for the leg area of the sling, increasing comfort and helping to prevent the material from roping and creasing under the resident's/patient's legs. Padded loop slings are available in four sizes (S, M, L and XL).

### **FEATURES**

- It provides good comfort and support for everyday residents/patients handling routines, such as bed and chair-related transfers.
- Padded sling with head support.
- Increased comfort and helps to prevent the material from roping and creasing under the resident's/patient's legs.
- Available in different sizes: S, M, L, XL

### **SPECIFICATIONS**

- Safe working load: 272 kg/600 lbs