



ARJO GENERAL PURPOSE POLYESTER SLING WITH HEAD SUPPORT



SKU: MLA2000-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Arjo standard padded loop sling is used on 2-point loop spreader bars. The slings are intended for everyday patient handling routines. The padded sling with head support offers padding for the leg area of the sling, increasing comfort and helping to prevent the material from roping and creasing under the resident's/patient's legs. Padded loop slings are available in four sizes (S, M, L and XL).

FEATURES

- It provides good comfort and support for everyday residents/patients handling routines, such as bed and chair-related transfers.
- Padded sling with head support.
- Increased comfort and helps to prevent the material from roping and creasing under the resident's/patient's legs.
- Available in different sizes: S, M, L, XL

SPECIFICATIONS

- Safe working load: 272 kg/600 lbs

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au