




## ARJO MAXI MOVE LIFTER



**SKU:** KMCSUN-D-22-H

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The Arjo Maxi Move Lifter is specifically engineered to facilitate the efficient and safe transfer and repositioning of patients or residents, even in situations that require significant physical effort. This innovative device empowers a single caregiver to handle the demanding tasks associated with daily patient mobility. To cater to a wide range of transfer requirements, the Arjo Maxi Move Lifter is equipped with a variety of spreader bars, allowing for easy customization and adaptation.

This flexibility ensures that the lift can accommodate the specific needs and conditions of each individual, providing a personalized and comfortable experience during transfers. By offering a comprehensive solution that enhances caregiver efficiency and patient comfort, the Arjo Maxi Move Lifter contributes to an improved standard of care and promotes the well-being of both caregivers and patients.

### KEY FEATURES

- **Stable Lifting:** Experience the utmost stability during lifts with the Arjo Maxi Move Lifter's consistent vertical lifting action, keeping the patient securely positioned at a constant distance from the mast.
- **Precision Transfers:** Optimize efficiency and minimize manual handling using the integrated powered dynamic positioning system (PDPS) of the Arjo Maxi Move Lifter, which seamlessly repositions the patient within the sling for precise transfers.
- **Versatility:** Customize the Arjo Maxi Move Lifter to individual needs by selecting from an extensive range of slings, made possible with the optional interchangeable spreader bars designed for adaptability.
- **Wide Lifting Range:** Effortlessly perform transfers onto beds or from the floor with the Arjo Maxi

Move Lifter's dual-action telescopic mast, which can be complemented by an extended jib for an expansive lifting range.

- Flexible System with Numerous Options Available: Enhance functionality and adaptability to different environments with optional features of the Arjo Maxi Move Lifter, such as an available scale for accurate weight measurements and three sizes of castors, including extra-low options for easy access under low couches and trolleys.

## SPECIFICATIONS

Weight and dimensions	Max safe working load
	kg
Max safe working load kg	227 kg
Max safe working load lb	500 lb
Max safe working load if equipped with extended jib kg	130 kg
Max safe working load if equipped with extended jib lb	287 lb
Lifter weight (with spreader bar and battery) kg	70 kg
Lifter weight (with spreader bar and battery) lb	154 lb
Max lifting height mm	1575 mm
Max lifting height in	62 in
Min lifting height mm	225 mm
Min lifting height in	8 ¾ in
Lifting stroke mm	1350 mm
Lifting stroke in	53,12 in
Electrical	
Lifter protection class	IPX4
Hand control protection class	IPX7
Battery	24V 4Ah

## BROCHURES

[Arjo Maxi Move Specifications Sheet](#)

[Maxi Move Instructions for use](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)