

ARJO MAXI TRANSFER SHEET

\$1,069.00



SKU: MAA7000R-GB-H



PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

Maxi Transfer Sheet is a repositioning and lifting sheet (which can also be used in place of bed linen) and is used with a passive lifting device for lateral and in-bed positioning. An exciting concept designed for lateral and inbed positioning solutions that simplify transfer routines.

FEATURES

STRAPS AND SEAMS

Colour-coded strap indicators at the edge of the sheet show the size of the sheet. All stitching is at the edge of the bed, ensuring a seam-free surface under the patient.

INTEGRATED SIDE VALANCES

The long valances cover the straps of the Maxi Transfer Sheet giving the appearance of a standard bed sheet.

SOFT FABRIC

Polyester fabric that is highly breathable and prevents heat build-up, to support microclimate management.¹

LOOP DESIGN

Soft loop fabric, and loop designed for attachment to the lifter of choice offering flexibility in repositioning the patient. Maxi Transfer Sheet is compatible with Maxi Move (floor lifter), as well as Maxi Sky 2 and Maxi Sky 2 Plus (ceiling lifters).

SPECIFICATIONS