


ARJO MAXI TRANSFER SHEET



SKU: MAA7000R-GB-H



PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

Maxi Transfer Sheet is a repositioning and lifting sheet (which can also be used in place of bed linen) and is used with a passive lifting device for lateral and in-bed positioning. An exciting concept designed for lateral and inbed positioning solutions that simplify transfer routines.

FEATURES

STRAPS AND SEAMS

Colour-coded strap indicators at the edge of the sheet show the size of the sheet. All stitching is at the edge of the bed, ensuring a seam-free surface under the patient.

INTEGRATED SIDE VALANCES

The long valances cover the straps of the Maxi Transfer Sheet giving the appearance of a standard bed sheet.

SOFT FABRIC

Polyester fabric that is highly breathable and prevents heat build-up, to support microclimate management.¹

LOOP DESIGN

Soft loop fabric, and loop designed for attachment to the lifter of choice offering flexibility in repositioning the patient. Maxi Transfer Sheet is compatible with Maxi Move (floor lifter), as well as Maxi Sky 2 and Maxi Sky 2 Plus (ceiling lifters).

SPECIFICATIONS

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriiffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au