



ARJO MAXI TWIN FLOOR LIFTER

\$9,809.00 - \$13,426.00 Price range: \$9,809.00 through \$13,426.00



SKU: KTBB4BSX2AU-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The twin mast design of the Arjo Maxi Twin Floor Lifter distributes the weight evenly, and the large handholds enable caregivers to manoeuvre the lifter. The open mast design also allows for additional legroom for the patient. The dual mast is designed to enable close, reassuring contact between caregiver and patient throughout the transfer. Integrated with the spreader bar is a powered dynamic positioning system (PDPS) that helps in repositioning the patient in the sling, minimising manual handling. Maxi Twin has been designed to be operated by a single caregiver.

Optional electronic scale enables the caregiver to weigh the patient while performing a transfer.

FEATURES

- Manoeuvrability
- The powered dynamic positioning system (PDPS) minimises manual handling and improves precision when positioning the patient
- Range of spreader bars
- Ergonomic handling
- The friction solution reduces swinging of the spreader bar during patient transfers
- Optional suspended scale supports an efficient weighing routine

SPECIFICATIONS

Maximum Lifting Height	1460mm
Minimum Lifting Height	230mm

Lifting Stroke	1230mm
Minimum Total Storage Height	1370mm
External Width Legs Closed	755mm
External Width Legs Open	1371mm
Internal Width Legs Closed	630mm
Internal Width Legs Open	1231mm
Leg Clearance Floor to Top of Legs Std Castors	115mm
Chassis Clearance Floor to Bottom of Chassis	31mm
Lifter Length	1298mm
Turning Diameter	1440mm
Weight with Spreader Bar and Battery	58kg
SWL	182kg