



ARJO MINSTREL PATIENT FLOOR LIFTER



SKU: HMA0007-AU-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Minstrel is the classic loop sling lifter for everyday transfer routines with dependent residents.

FEATURES

- Easy to understand and simple to operate
- The handset enables the caregiver to move and stay close to the resident and maintain good contact throughout the lift and transfer process
- A single caregiver can move a resident who has fallen on the floor
- Available with an optional scale that enables the caregiver to weigh the resident while performing a transfer

SPECIFICATIONS

| Maximum Lifting Height | 2000mm |
|--------------------------------------|--------|
| Minimum Lifting Height | 625mm |
| Lifting Stroke | 1375mm |
| Minimum Total Storage Height | 1330mm |
| External Width Legs Closed | 645mm |
| External Width Legs Open | 1181mm |
| Internal Width Legs Closed | 545mm |
| Internal Width Legs Open | 1080mm |
| Leg Clearance Floor to Top of Legs | 110mm |



| Chassis Clearance Floor to Bottom of Chassis | 45mm |
|--|--------|
| Lifter Length | 1275mm |
| Turning Diameter | 1420mm |
| Weight with Spreader Bar and Battery | 56kg |
| SWL | 190kg |

For further information go to medimart.com.au or give us a call: 1300 633 063

SA

Brighton Showroom

343 Brighton Rd, North Brighton SA 5048 08 7079 8773 brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112 08 7079 8777 elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd, Reynella SA 5161 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199 03 8763 3451 frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168 03 7503 8354 monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820 1300 633 468 support@medimarnt.com.au