



## ARJO MINSTREL PATIENT FLOOR LIFTER

**\$5,793.00 - \$9,136.00** Price range: \$5,793.00 through \$9,136.00



**SKU:** HMA0007-AU-H

## PRODUCT DESCRIPTION

## PURPOSE & BENEFITS

Minstrel is the classic loop sling lifter for everyday transfer routines with dependent residents.

## FEATURES

- Easy to understand and simple to operate
- The handset enables the caregiver to move and stay close to the resident and maintain good contact throughout the lift and transfer process
- A single caregiver can move a resident who has fallen on the floor
- Available with an optional scale that enables the caregiver to weigh the resident while performing a transfer

## SPECIFICATIONS

<b>Maximum Lifting Height</b>	2000mm
<b>Minimum Lifting Height</b>	625mm
<b>Lifting Stroke</b>	1375mm
<b>Minimum Total Storage Height</b>	1330mm
<b>External Width   Legs Closed</b>	645mm
<b>External Width   Legs Open</b>	1181mm
<b>Internal Width   Legs Closed</b>	545mm
<b>Internal Width   Legs Open</b>	1080mm
<b>Leg Clearance   Floor to Top of Legs</b>	110mm

<b>Chassis Clearance   Floor to Bottom of Chassis</b>	45mm
<b>Lifter Length</b>	1275mm
<b>Turning Diameter</b>	1420mm
<b>Weight with Spreader Bar and Battery</b>	56kg
<b>SWL</b>	190kg